

<b>Table of contents</b>
--------------------------

\\MedPhys
-----------

<table border="1"><tr><td>Technologies</td></tr></table>	Technologies
Technologies	

<table border="1"><tr><td><table border="1"><tr><td>Qingping</td></tr></table></td></tr></table>	<table border="1"><tr><td>Qingping</td></tr></table>	Qingping
<table border="1"><tr><td>Qingping</td></tr></table>	Qingping	
Qingping		

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td>dataacq_imagerecon_workflow</td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td>dataacq_imagerecon_workflow</td></tr></table></td></tr></table>	<table border="1"><tr><td>dataacq_imagerecon_workflow</td></tr></table>	dataacq_imagerecon_workflow
<table border="1"><tr><td><table border="1"><tr><td>dataacq_imagerecon_workflow</td></tr></table></td></tr></table>	<table border="1"><tr><td>dataacq_imagerecon_workflow</td></tr></table>	dataacq_imagerecon_workflow	
<table border="1"><tr><td>dataacq_imagerecon_workflow</td></tr></table>	dataacq_imagerecon_workflow		
dataacq_imagerecon_workflow			

<table border="1"><tr><td><table border="1"><tr><td><table border="1"><tr><td><a href="#">product_t1_mprage_sag_p2</a></td></tr></table></td></tr></table></td></tr></table>	<table border="1"><tr><td><table border="1"><tr><td><a href="#">product_t1_mprage_sag_p2</a></td></tr></table></td></tr></table>	<table border="1"><tr><td><a href="#">product_t1_mprage_sag_p2</a></td></tr></table>	<a href="#">product_t1_mprage_sag_p2</a>
<table border="1"><tr><td><table border="1"><tr><td><a href="#">product_t1_mprage_sag_p2</a></td></tr></table></td></tr></table>	<table border="1"><tr><td><a href="#">product_t1_mprage_sag_p2</a></td></tr></table>	<a href="#">product_t1_mprage_sag_p2</a>	
<table border="1"><tr><td><a href="#">product_t1_mprage_sag_p2</a></td></tr></table>	<a href="#">product_t1_mprage_sag_p2</a>		
<a href="#">product_t1_mprage_sag_p2</a>			

\\MedPhys\Technologies\Qingping\dataacq\_imagerecon\_workflow\product\_t1\_mprage\_sag\_p2

TA: 5:38 min Coil Selection: Manual Voxel Size: 1.0×1.0×1.0 mm<sup>3</sup> Acc:: 2 Rel. SNR: 1.00**Properties**

Start measurement without further preparation	On
Wait for User to Start	Off
Start measurements	Single Measurement
Prio Recon	Off
Auto Open Inline Display	Off
Auto Close Inline Display	Off
Load Images to MR View&GO	On
Auto Store Images	On
Disable auto transfer to PACS	Off
Load Images to Stamp Segments	Off
Load Images to Graphic Segments	Off
Graphic segment	Default
Inline Movie	Off

**Routine**

Slab Group	1
Slabs	1
Distance Factor	50 %
Position	L0.0 P14.4 F2.2 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	192
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	256 mm
FOV Phase	93.8 %
Slice Thickness	1.00 mm
TR	2500.0 ms
TE	3.37 ms
Averages	1
Concatenations	1
AutoAlign	Head > Basis
Coil Elements	HE1-4

**Contrast - Common**

TR	2500.0 ms
TE	3.37 ms
Magn. Preparation	Non-sel. IR
TI	1100 ms
Flip Angle	7 deg
Fat-Water Contrast	Standard
Dark Blood	Off
Contrasts	1
Reconstruction	Magnitude

**Contrast - Dynamic**

Dynamic Mode	Standard
Measurements	1
Multiple Series	Off
Reordering	Linear

**Resolution - Common**

FOV Read	256 mm
FOV Phase	93.8 %
Slice Thickness	1.00 mm
Base Resolution	256
Phase Resolution	100 %
Slice Resolution	100 %
Interpolation	Off

**Resolution - Acceleration**

Acceleration Mode	GRAPPA
Reference Scans	Integrated
Acceleration Factor PE	2
Reference Lines PE	32
Acceleration Factor 3D	1
Phase Partial Fourier	Off
Slice Partial Fourier	Off
Asymmetric Echo	Off
Elliptical Scanning	Off

**Resolution - Filter**

Raw Filter	Off
Elliptical Filter	Off
Distortion Correction	3D
Normalize	Off
Image Filter	Off

**Geometry - Common**

Slab Group	1
Slabs	1
Distance Factor	50 %
Position	L0.0 P14.4 F2.2 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
Slices per Slab	192
Phase Oversampling	0 %
Slice Oversampling	0.0 %
FOV Read	256 mm
FOV Phase	93.8 %
Slice Thickness	1.00 mm
TR	2500.0 ms
Multi-Slice Mode	Single Shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab Group	1
Position	L0.0 P14.4 F2.2 mm
Orientation	Sagittal
Phase Encoding Dir.	A >> P
AutoAlign	Head > Basis

**Geometry - AutoAlign**

Initial Position	L0.0 P14.4 F2.2
L	0.0 mm
P	14.4 mm
F	2.2 mm
Initial Orientation	Sagittal
Initial Rotation	0.00 deg

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table Position	2 mm
Table Position	F
Inline Composing	Off

**System - Miscellaneous**

Coil Selection	Manual
Radial Sorting	Off
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combination	Sum of Squares
Matrix Optimization	Off

**System - Adjustments**

Adjustment Strategy	Standard
B0 Shim	Standard
B1 Shim	TrueForm
Adjustment Tolerance	Auto
Adjust with Body Coil	Off
Confirm Frequency	Never
Assume Silicone	Off

**System - Adjust Volume**

Position	L0.0 P14.4 F2.2 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	240 mm
F >> H	256 mm
R >> L	192 mm
Reset	Off

**System - pTx**

B1 Shim	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.194268 MHz
? Ref. Amplitude 1H	0.000 V
Reset	Off
Correction Factor	1.00
Image Scaling	0.800

**Physio - Signal**

1st Signal/Mode	None
TR	2500.0 ms
Concatenations	1

**Physio - Cardiac**

Fat-Water Contrast	Standard
Magn. Preparation	Non-sel. IR
TI	1100 ms
Dark Blood	Off
FOV Read	256 mm
FOV Phase	93.8 %
Phase Resolution	100 %

**Physio - PACE**

Resp. Control	Off
Concatenations	1

**Inline - Subtraction**

Subtract	Off
Measurements	1
StdDev	Off
Save Original Images	On

**Inline - MIP**

MIP Sag	Off
MIP Cor	Off
MIP Tra	Off
MIP Time	Off
Radial MIP	Off
Save Original Images	On
MPR Sag	Off
MPR Cor	Off
MPR Tra	Off

**Inline - Composing**

Inline Composing	Off
------------------	-----

**Inline - MapIt**

MapIt	None
Flip Angle	7 deg
Measurements	1
Contrasts	1
TE	3.37 ms
TR	2500.0 ms
Save Original Images	On

**Sequence - Part 1**

Sequence Name	tfl
Dimension	3D
Excitation	Non-sel.
RF Pulse Type	Fast
Gradient Mode	Normal
Flow Compensation	None
Reordering	Linear

**Sequence - Part 1**

Bandwidth	200 Hz/Px
Echo Spacing	7.94 ms
Asymmetric Echo	Off
Turbo Factor	192

**Sequence - Part 2**

Introduction	Off
RF Spoiling	On
Incr. Gradient Spoiling	Off

**Sequence - Assistant**

SAR Assistant	Off
---------------	-----